

TK-VETKIT

CE


VEC-15 Series

Portable Veterinary Endoscopes



Clinical Applications

- ▶ VEC-15 series are suitable for animals from small to big, enabling precise examination and treatment of key areas like the stomach, colon, esophagus, uterus, respiratory tract, and ENT. This series is adaptable to various applications especially for scenarios requiring portability.



Product Features



Package Solution



Image Processor (VHU-1405)

- 4.0-inch touch screen
- 64G TF Card
- Rechargeable battery supporting 2.5 working hours

Endoscope (VEC-15 Series)

- High-Definition endoscope
- Convenient for use anytime, anywhere

Portable Irrigation and Insufflation System (VOTH-003)

- Adjustable maximum air pump pressure range: 35-55kPa

Specifications

Model	Insertion Tube	Working Channel	Working Length	Field of View	Depth of Field	Bending (Up/Down)	Rotation (Left / Right)	Function Keys	Water / air nozzle	Water Jet
VEC-15	9.2mm	3.2mm	1500mm	120°	3-100mm	180°/130°	160°/160°	3	Yes	\
VEC-15L	9.2mm	3.2mm	3000mm	120°	3-100mm	180°/130°	160°/160°	3	Yes	\

Add: Room 301, 3F, Plant 1, No. 199, Shui'an 1st Road, Xiangzhou District, Zhuhai, China.
Tel: +86-755-23761871 E-mail: intl@seesheen.com Web: www.seesheen.com

YouTube TKVETKITENDOSCOPE

in TK VETKIT

Instagram tkvetkit520

ISO 13485